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Non-contact Infrared Thermometer
User Manual
Model No. BM-041

Product Description

This product is a non-contact multi-functional infrared thermometer, could have a professional measurement of human temperature, using a German-made high-performance temperature sensor, can accurately measure body temperature in different environment, the notice of LED display and sound is very easily for users to read degrees and judge the measurement results. There is a LED display, can simply and easily to read the degrees, measuring your body temperature in just one second, so it is a high-quality product you can trust.

Scope of application

The body temperature or surface temperature of the object will be displayed through non-contact measuring thermal radiation on the surface of a temple or object.

Basic Principle

When all objects have a temperature higher than absolute zero, they radiate a certain amount of infrared energy, different categories of objects radiated out of the infrared energy value is different, according to the radiation of infrared energy and wavelength can determine the surface temperature value of such objects. According to this principle, the temperature of the temple position on both sides of the forehead can be measured accurately because of its relative stability. At the same temperature, the infrared intensity emitted by the temple is not the same as the infrared energy intensity radiated by other objects (e.g. water), so there are at least two measurement modes (body temperature mode and object temperature mode). The Object temperature mode is the measurement the temperature of pure water as a benchmark, so when you measure other objects, there may be some deviation due to the different properties of the object, which is normal.

Product Functions

- Non-contact measurement
2~3cm distance from the surface of the temple
Avoid cross infection
- German-made sensor
The precision of the body temperature is ≤±0.2℃
- LED alarm for body temperature
Easily to measure in the dark environment
Body temperature mode (≤37.9℃ LED white color, ≥38.0℃ LED Red color)
- Fever alarm
Major for baby and children
Above or equal to 38℃ automatic sound and LED red alarm
- Multi-function
Body temperature, milk temperature, bath water temperature, room temperature
Take care of the baby's healthy life
- Body/Object temperature (milk, bath water, room) 2 modes
Celsius/Fahrenheit units can choose
30 seconds no operation automatic shut down


Notice before use

In order to obtain stable and reliable measurement data, when the environment temperature changes significantly, the product should be placed indoors after 15-20 minutes before use. Depending on the skin color, skin thickness and body parts, the temperature measured will be different, which is normal. The more exposed the body parts of the body will be the greater the influence of environment temperature. Also don't expose the product in the sun and don't touch the water.

Reference range

Body parts	Normal body temperature	Normal body temperature
Anus	36.6-38.0	97.8-100.4
Oral	35.5-37.5	95.9-99.5
Armpit	35.5-37.3	95.9-99.1
Ear	35.8-38.0	96.4-100.4
Temple	36.0-37.9	96.8-100.2

Normal body temperature range at different ages


* Children's body temperature is higher

Specifications


Type	Non-contact
Effective distance	2~3cm from the temple
Measuring range	Body 34.0℃~42.9℃ (93.2℉~109.2℉)
	Object 10.0℃~70.0℃ (50.0℉~158.0℉)
Precision	Body ±0.2℃ (±0.4℉)
	Object ±1.0℃ (±1.8℉)
Resolution	0.1℃/℉
Operate condition	16℃-35℃ (60.8℉-95.0℉) RH≤85%
Storage condition	-20℃-50℃ (-4.0℉-122.0℉) RH≤85%
Power	DC3V, AA Battery x 2
Temperature Unit	Celsius ℃ and Fahrenheit ℉F
Auto power off	Within 30 seconds
Net Weight	75g (included battery)

Unit Parts Introduction

The non-contact multi-function infrared thermometer consists of an infrared temperature sensor, LED display, battery chamber, keystrokes, mode switch and main case, etc., as shown below:

Keys Operation


Displays 👁 means body temperature mode
Displays 🏠 means object temperature mode (room/bath/water/milk)

CF/℉ switch

How to use

Initial use or just after you put in a new battery, preheat for about ten minutes to make the measurement data more accurate. Switch the product to body temperature mode (when the human head symbol displayed before starting to measure), make sure that the product is in body temperature mode, maintain a distance of 2-3cm perpendicular to the skin of temple, then press the measuring key and hear a beep deep to indicate that the measurement is complete, at this moment you can remove the product from the temple site and can read the measured value on the LED screen. Do not be hurried and wait for the measurement to be completed, otherwise it will be incorrect measurements. If the product is not operated within 30 seconds after measurement, it will be automatically shutdown.

1. Measure body temperature

Measure temple if no sweating
Measure back of ear if sweating

2. Measure milk, bath water or room temperature

Milk
Room
Bath Water

Warning

- Please read this use instruction manual carefully before use.
- The environment temperature to use this product is 16℃ to 35℃, the optimum temperature is 25℃.
- The accuracy of this product is greatly reduced when used in environments where the ambient temperature is above 35℃ or below 16℃.
- Do not place the product near charged objects to avoid electric shock.
- Do not use this product in an environment where the relative humidity is higher than 85%.
- Do not bring the product too close to the electromagnetic field (e.g. radios, mobile phones, etc.).
- Do not place the product to expose in the sun or near the stove, and stay away from water.
- This product is a precision device, do not disassemble, twist, bite, fall, impact.
- If the hair, sweat, hat or scarf cover the temple, that will affect the accuracy of the measurement results, by the way, please confirm that the measuring distance is 2 to 3 cm from the skin of temple.
- In order to obtain accurate, stable and reliable measurements, use the product after 10 minutes indoors when room temperature varies significantly.
- There are significant changes in the environment before and after the two measurements (e.g. from outdoor to indoor), please wait 30 minutes in the new environment before measuring.
- It is dangerous for patients to judge and treat their condition solely on the basis of temperature measurements. Be sure to follow your doctor's instructions.
- Do not throw the product and the attached battery into the fire and burn it.
- In the event of a problem with this product, contact local distributor and do not attempt to repair it yourself.
- This product is not a toy, please place it out of the reach of children.

Warranty Certificate

Please fill out this warranty certificate and fax /e-mail with the invoice's copy within 20 days after purchase. Please take a copy of the invoice if it is a thermal paper.

We guarantee to offer one year free service from the date of purchase within the limit of the following conditions:
 (1) For any defect which in our judgment is caused under normal use, we are responsible for repairing or replacing parts excluding cabinet, consumptive parts and accessories. Free of charge within the one year guarantee period commencing from the date of purchase. Any defective part which has been replaced shall become our property.
 (2) This guarantee certificate shall automatically be void if (a) it has been explicitly or implicitly modified, tempered with altered or repaired in any way by persons not authorized by our company, or (b) it has been damaged through misuse, negligence, liquid ingress, accident, insects or natural calamities, or (c) it has become manufactured due to power interruption or abnormal voltage input.
 (3) For service, the user will be responsible for the transportation of the product to and from our service stations for repair.
 (4) This guarantee certificate is valid only in Hong Kong and Macau.
 We will charge for the parts and workman fee after one year guarantee.
 In case there is any inconsistency or conflict between the English version and the Chinese version of the terms, the Chinese version shall prevail.

Name: _____
 Product: Babymate non-contact thermometer
 Model no.: BM-041
 Purchase date: _____
 Purchase shop: _____
 email address: _____

Golden Way International Enterprises Ltd.
 Address: Flat 2, 15/F, Goldfield Ind. Cr., 1st Wai Wo Rd., Kowloon, N.T., Hong Kong
 Tel: 2666-2760 Fax: 2666-3370
 Email: service@golden-way.com.hk
 Web: www.golden-way.com.hk